


DISPOSABLE BOUGIE DILATOR SETS



Indicated for Oesophageal Dilatation

Features:

- Two radiopaque markers on each end for increased visibility under x-ray
- Markings on both sides for easy identification of insertion depth
- Conical angle at tip allows for ease of insertion
- Made from medical grade polymer

Specyfication:			
Model	Dilator Diameter(mm)	Working Length(mm)	Guide Wire
BDS - 05	5	850	0.035"
BDS - 06	6	850	0.035"
BDS - 07	7	850	0.035"
BDS - 08	8	850	0.035"
BDS - 09	9	850	0.035"
BDS - 10	10	850	0.035"
BDS - 11	11	850	0.035"
BDS - 12	12	850	0.035"
BDS - 13	13	850	0.035"
BDS - 14	14	850	0.035"
BDS - 15	15	850	0.035"
BDS - 16	16	850	0.035"
BDS - 17	17	850	0.035"
BDS - 18	18	850	0.035"
BDS - 19	19	850	0.035"
BDS - 20	20	850	0.035"

