


UROFINO

SINGLE USE CYSTOSCOPE

- ▶ Tapered-Shape distal end aids insertion
- ▶ Dual light guide enhances visualisation
- ▶ Larger working channel available for added flexibility
- ▶ HD image quality on par with conventional scopes



Universal Processor - compatible with range of disposable scopes



Available in 2 sizes with the option of larger working Channel

Model	UV-100-E	UV-100-T-E
Field of View	110°	110°
Depth of View	5-50mm	5-50mm
Deflection Up/Down	210°/130°	210°/130°
Insertion Diameter	4.7mm	5.3mm
Distal End Diameter	4.7mm	5.3mm
Working Length	38.5cm	38.5cm
Working Channel	2.2mm	2.6mm

Brought to you by:

CE 0123 EN ISO13485:2016

SEEGEN 视介光电
vision borderless

Benefits of Urofino

Single-Use Cystoscope

1

Safety - A new sterile scope for every patient, reducing the risk of infection

2

Consistency - Predicable clinical performance everytime - no issue with wear and tear

3

24 x 7 Availability - No need to wait for reprocessing of reusable scopes

4

Cost Transparency - No Service or Leasing Agreements

5

Maintenance Free - **No** degradation, **No** reprocessing and **No** repairs

6

Compatibility - works with most commonly used diagnostic and therapeutic tools

Distributed in Ireland by:

FLEETWOOD
HEALTHCARE
EXCELLENCE THROUGH
INNOVATION



Fleetwood Healthcare

ROI: Unit 5 Bromley Business Park, Kilpedder, Co. Wicklow, A63X924
Tel: +353 1 201 4620

NI: Innis House, Holywood, Co. Down, BT18 9HF
Tel: +44 28 9099 5160

e-mail: sales@fleetwoodhealthcare.com

Website: www.fleetwoodhealthcare.com